



ERYTHROCYTE SEDIMENTATION RATE ANALYZERS

The Most Efficient ESR Solutions for Your Laboratory



TRANSFORMING ONE OF THE WORLD'S MOST COMMONLY ORDERED LAB TESTS


Sed rates are one of the oldest and most common tests in the world. Isn't it time they got easier to perform? ALCOR® Scientific's iSED® analyzers are designed to make ESR testing as easy and efficient as possible, providing results within 20 seconds with the least amount of hands-on time.



Results within 20 seconds



Walk-away automation



Highly correlated to the Westergren Method



Sampling from capped primary EDTA tube



100 µl aspirated volume


ESR RESULTS IN SECONDS

ALCOR Scientific's iSED technology utilizes photometric rheology to measure the rouleaux formation, the earliest and most critical phase of red blood cell sedimentation². An EDTA blood sample is injected into the iSED analyzer's flow cell which causes the red blood cells to disaggregate. The flow cell sensors capture the difference in light transmission as the red blood cells begin to re-aggregate.

The technical innovation of iSED ESR analyzers consists of "directly" measuring the aggregation of the red blood cells, whereas the traditional methods "indirectly" measure the aggregation of the red blood cells by recording the length at which the red blood cells settle in a Westergren tube.

The iSED technology is highly correlated to the Westergren method and provides results in mm/hr. Test results are available in just 15-20 seconds without the drawbacks and inconvenience of traditional ESR testing.

Flow Cell:




1. ALCOR Scientific internal data.

2. Fabry, T. L. Mechanism of erythrocyte aggregation and sedimentation. *Blood*. 1987;70(5):1572-1576.



SCARCE LABORATORY RESOURCES DEMAND THE MOST EFFICIENT ESR WORKFLOW


SINGLE TUBE WORKFLOW



Primary EDTA sample tube



2




< 3 seconds
of hands on time
per sample¹

iSED ELITE ESR Analyzer




3



ESR Results

SAMPLE RACK WORKFLOW




Primary EDTA sample tubes
in a hematology rack



< 1 second
of hands on time
per sample¹

2



iSED PRO ESR Analyzer



3



ESR Results

Results within
20 seconds

QUALITY CONTROL MADE SIMPLE

SEDiTROL® Quality Controls are bi-level, human-based whole blood controls in ready-to-use, barcoded tubes. All SEDIROL customers have access to iQAP, ALCOR Scientific's online peer-to-peer quality assurance program. No special preparation or washing is required to run quality control.



WHICH ANALYZER IS
RIGHT FOR YOU?



iSED[®] PRO

Best for higher volume labs that prefer a fully automated, rack-based workflow.



iSED[®] ELITE

Best for medium-high volume labs that prefer a fully-automated, flexible, single tube workflow.



mini iSED[®]

Best for lower volume labs that prefer a compact analyzer footprint.

Sample Capacity	120 samples (12 Sysmex Hematology Racks)	20 samples	1 sample
Test Throughput	Up to 180 results/hour (Rack based)	Up to 180 results/hour	Up to 180 results/hour
Sample Volume	500 µl total volume (100 µl aspirated volume)	500 µl total volume (100 µl aspirated volume)	500 µl total volume (100 µl aspirated volume)
Sample Type	Primary EDTA tube; pediatric tube compatible, samples must be in Sysmex hematology racks	Primary EDTA tube; pediatric tube compatible	Primary EDTA tube; pediatric tube compatible
Test Time	Results within 20 seconds	Results within 20 seconds	Results in 15 seconds
Quality Control	Onboard storage and automatic scheduling Room temperature storage, 60 day open vial stability	User initiated control runs Room temperature storage, 60 day open vial stability	User initiated control runs Room temperature storage, 60 day open vial stability
Barcode Reader	Internal	Internal	Internal
LIS Connectivity	Serial RS232 DB9 port; bi-directional	LAN TCP/IP or Serial RS232 DB9 port; uni-directional	LAN/IP uni-directional
Dimensions (L x W x H)	49.5 x 55.3 x 67 (cm) 20 x 22 x 27 (in)	36 x 27 x 35 (cm) 14.17 x 10.63 x 13.78 (in)	36 x 19 x 24 (cm) 14 x 7.5 x 9.5 (in)
Weight	38 kgs; 84 lbs	13.6 kgs; 30 lbs	4.5 kgs; 10 lbs
Printer	External (Optional)	Internal	External (Optional)



Product Information		
Catalog #	Description	Unit
112-00120-SYS	iSED PRO ESR Analyzer Series S	1 each
112-00222	iSED ELITE ESR Analyzer	1 each
112-01017	miniISED® ESR Analyzer	1 each
112-00250	Test Card; 250 test credits	1 each
112-01000	Test Card; 1,000 test credits	1 each
112-02000	Test Card; 2,000 test credits	1 each
112-05000	Test Card; 5,000 test credits	1 each
120-07500	iSED PRO Test Card; 7,500 test credits	1 each
120-12500	iSED PRO Test Card; 12,500 test credits	1 each
120-25000	iSED PRO Test Card; 25,000 test credits	1 each
120-50000	iSED PRO Test Card; 50,000 test credits	1 each
DSC06	SEDITROL Quality Control (Levels 1 & 2)	Pack of 6 x 4.0 mL
DSC01		Pack of 2 x 4.0 mL
112-12-020	deepCLEAN® Cleaning Solution for miniISED, iSED ELITE, and iSED	3 x 4.0 mL
112-12-022	deepCLEAN PRO Cleaning Solution for iSED PRO	3 x 4.0 mL
112-12-006	iWASH® PRO Wash Fluid	1 x 1.8 L
112-12-001	iWASH Wash Fluid	Pack of 4
112-12-003	miniWASH® Wash Fluid	Pack of 4
112-12-009	iWASTE® PRO Container	1 x 1.8 L
112-12-002	iWASTE Container	Pack of 24
112-12-005		Pack of 4
112-12-004	miniWASTE® Container	Pack of 4
DS-05233	ALCOR Thermal Printer Paper	Pack of 5
DS-05240	ALCOR Thermal Printer	1 each



miniISED®

iSED® ELITE

The most
versatile family
of ESR solutions!

iSED® PRO

ORDERING

ALCOR Scientific Sales: 1.800.495.5270
info@alcorscientific.com

SERVICE AND SUPPORT

ALCOR Scientific Technical Support: 1.800.495.5270
techservice@alcorscientific.com

20 Thurber Boulevard
Smithfield, RI, 02917
USA
401.737.3774
ALCORSCIENTIFIC.COM



ALCOR, iSED, miniISED, SEDITROL, deepCLEAN, iWASH, iWASTE, miniWASH, and miniWASTE are registered trademarks of ALCOR Scientific LLC.