



URGENT CARE

Expand your test menu with 20 second sed rates.

WHY NOT OFFER SED RATE TESTING?



Expand your test offering and generate revenue with an easy-to-use screening test



Improve patient satisfaction with a test that can be performed during a patient visit



Improve workflow synergies by doing more with a lavender top tube



INTRODUCING RAPID, COMPACT, & EFFICIENT ESR ANALYZERS FROM ALCOR SCIENTIFIC

The miniSED and iSED analyzers generate an ESR result within 20 seconds from a capped EDTA sample tube. Each ESR test only requires 100 μ l of sample to be aspirated and approximately 3 seconds of hands-on time.



PRODUCT DETAILS

iSED[®] ELITE



Sample Capacity	20 samples (random access loading)
Sample Type	Primary EDTA tube; pediatric tube compatible
Sample Volume	500 µl total volume (100 µl aspirated volume)
Test Throughput	Up to 180 results / hour
Sample Mixing	Automatic onboard
Test Time	Results within 20 Seconds
Barcode Reader	Internal
Dimensions (L x W x H)	36 x 27 x 35 (cm) 14.17 x 10.63 x 13.78 (in)
Data Storage Capacity	>4 million
LIS Connectivity	LAN TCP/IP or Serial RS232 DB9 port
Printer	Internal

miniSED[®]



Sample Capacity	1
Sample Type	Primary EDTA tube; pediatric tube compatible
Sample Volume	500 µl total volume (100 µl aspirated volume)
Test Throughput	Up to 180 results / hour
Sample Mixing	External before sample is loaded
Test Time	Results in 15 seconds
Barcode Reader	Internal
Dimensions (L x W x H)	36 x 19 x 24 (cm) 14 x 7.5 x 9.5 (in)
Data Storage Capacity	>4 million
LIS Connectivity	LAN TCP/IP
Printer	Optional external printer

ORDERING INFORMATION

Catalog #	Description	Unit
112-00222	iSED ELITE ESR Analyzer	1 each
112-01017	miniSED ESR Analyzer	1 each
112-00250	Test card; 250 test credits	1 each
112-01000	Test card; 1,000 test credits	1 each
112-02000	Test card; 2,000 test credits	1 each
112-05000	Test card; 5,000 test credits	1 each
DSC06	SEDiTROL Quality Control (Levels 1 & 2)	Pack of 6 x 4.0 mL
DSC01		Pack of 2 x 4.0 mL
112-12-001	iWASH [®] Wash Fluid	Pack of 4
112-12-003	miniWASH [®] Wash Fluid	Pack of 4
112-12-002	iWASTE [®] Container	Pack of 24
112-12-005		Pack of 4
112-12-004	miniWASTE [®] Container	Pack of 4
DS-05233	ALCOR Thermal Printer Paper	Pack of 5
DS-05240	ALCOR Thermal Printer	1 each

20 Thurber Boulevard
Smithfield, RI, 02917
USA
401.737.3774
ALCORSCIENTIFIC.COM

ORDERING

ALCOR Scientific sales: 1.800.495.5270

info@alcorscientific.com

ALCOR, iSED, and miniSED are registered trademarks of ALCOR Scientific LLC.

1017-09-201, Rev. 0

