# SENTINEL plus®

Enteral Feeding Pump

Become part of the evolution



Delivering the same durability & reliability, simply smaller.

Now that's a plus!







#### **DUAL DISPLAY**

The bright LED display allows caregivers to simultaneously and conveniently view the infusion rate and volume delivered from a distance.



### **ADVANCED MEMORY**

The pump operating history is retained until manually erased, even in the event of power failures or when pump is turned off.



#### **DOSE CHECK**

A record is kept of the accumulated amount of formula delivered over multiple feedings.



#### SIMPLE ALARMS

Audible and visual alarms are easily differentiated.



#### **CLOVER-SHAPED ROTOR**

The unique clover-shaped rotor allows for easy recognition of pump operation from a reasonable distance.



## SENTINELplus® Advanced Technology

The SENTINELplus® Enteral Feeding Pump is a revised version of our popular SENTINEL® product. The lighter, more compact design brings the same simplicity and reliability of its predecessor. Contact us to learn more!

#### **Order Information**

Cat.

Description

S-3000-IU

SENTINELplus®, **Enteral Feeding Pump**  **Packaging** 

1 each

**Specifications** 

Mechanism Rotary Peristaltic

5-295 mL/hour in 1 mL increments **Operating Range** 1-2000 mL in 1 mL increments **Dose Limit** 

2000-9999 in 5 mL increments

± 10% of selected flow rate with ALCOR sets Dose Accuracy Occlusion Pressure

Battery

Does not exceed 15 psi

Battery backup 6 hours minimum at 125

mL/hr

**Dimensions** 15cm H x 20cm W x 9cm D

Weight 1.42kg

**Power Options** All country-specific connectors available

> www.alcorscientific.com info@alcorscientific.com 401.737.3774 800.495.5270

> > 113-09-500 11/2016