

# DELIVERING THE SAME DURABILITY & RELIABILITY, SIMPLY SMALLER. NOW THAT'S A PLUS!







# **DUAL DISPLAY**

The bright LED display allows caregivers to simultaneously and conveniently view the infusion rate and volume delivered from a distance.



# **ADVANCED MEMORY**

The pump operating history is retained until manually erased, even in the event of power failures or when pump is turned off.



### **DOSE CHECK**

A record is kept of the accumulated amount of formula delivered over multiple feedings.



# **SIMPLE ALARMS**

Audible and visual alarms are easily differentiated.



### **CLOVER-SHAPED ROTOR**

The unique clover-shaped rotor allows for easy recognition of pump operation from a reasonable distance.

### **SPECIFICATIONS**

Mechanism Rotary Peristaltic

Operating Range 5-295 mL/hour in 1 mL increments

Dose Limit 1-2000 mL in mL increments

2000-9999 in 5 mL increments

Dose Accuracy + 10% of selected flow rate with ALCOR sets

Occlusion Pressure Does not exceed 15 psi

Battery Battery backup 6 hours minimum at 125 ML/hour

**Dimensions** 15cm H x 20cm W x 9cm D

Weight 1.42kg

Power Options All country-specific connectors available





